# CareIP Mobile

Bruksanvisning | User Guide | Benutzerhandbuch

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## Read first: Safety information

- The compartment covers on the reverse may be opened only by authorized persons in accordance with the instructions in the service manual.
- Battery replacement may be performed only by authorized persons in accordance with the instructions in service manual. Only use recommended battery type as stated in the service manual. CAUTION – Risk of explosion if battery is replaced by an incorrect type. Dispose of used batteries according to the instructions.
- The antenna shall always be mounted if the unit is in operation.
- Installation must be performed only by authorized persons in accordance with the instructions in the service manual. Only use recommended power supply as stated in the service manual.

Always read and follow the safety information accompanied by this symbol

## Important Information

All systems using radio and network communication are subject to interference beyond the user's control.

Products from CareTech are designed to minimise the impact of such interference. Nevertheless, the user must be aware that system components can be subjected to interference or other influences that may cause malfunction.

It is therefore important to regularly check that every part of the system works in all areas, especially radio communications. Contact your supplier immediately in case of any suspected malfunction.

Users should pay particular attention to the risk of disruption from products which communicate using the same or adjacent frequencies.


For further information, please contact your supplier or refer to the service manual or the documentation at [www.caretech.eu](http://www.caretech.eu)

## Package contents

- CareIP Mobile carephone
- Elliot radio trigger with accessories
- Network cable (cat5E UTP TIA568B)
- Power supply
- Pull cord
- User guide

## CareIP parts

There are three buttons on the front, with the following standard functions:

The following indicator lights can be found on the front of the CareIP Mobile:

### On, (Green)

- Steady light: Mains power.
- Flashing light: Mains failure or radio programming mode.

### Alarm, (Red)

- Steady light: Idle mode.
- Slow flashing: Alarm state.


### Error, (Red)

- Off: No faults detected.
- Steady light: Problem with Ethernet link or IP configuration.  
(NOTE! Always steady light if only GSM connection is used).
- Slow flashing: Registration to SIP server failed.

### Status, (Red/Green)

- Off: No faults detected.
- Steady red light: No connection to GSM network.

The following can be found on the reverse of the CareIP Mobile.

 Volume control: Loudspeaker volume adjustable in eight levels.

**PC** Pull cord: Attachment for pull cord when wall mounted.

 Power switch: Turns CareIP ON and OFF.

 DC IN: Connection for power supply 7,5VDC

 Ethernet port: Connection for Ethernet cable from broadband/network.

 I/O connectors (wired). Wiring must be performed only by authorized persons in accordance with the instructions in the service manual.

With CareIP GSM a standard antenna is delivered, if other antenna will be used it must be adapted for use with CareIP Mobile.

 The compartment covers on the reverse may be opened only by authorized persons in accordance with the instructions in the service manual.

## **Description of functions**

CareIP Mobile is a digital alarm unit with a built in GSM module, the unit has possibility to transfer alarms over broadband/network connections and mobile telephone networks. The unit offers the user the possibility of two-way voice communication using existing standard for VoIP (Voice over Internet Protocol) calling. The unit uses SIP, a standard protocol for transfer of speech and video over broadband connections. In addition to voice-connected alarm calls CareIP offers a platform for transfer of detailed alarm information, speech messages, medical information, online status etc.

The Elliot radio trigger is adapted for use with the unit. Elliot can be worn using a wristband, neckband or pocket clip. The CareTech Easy Press Concept (EPC) accessory makes it easier for those with reduced physical strength to activate the alarm.

### **Activating the alarm**

The assistance alarm can be activated with a radio trigger or by pressing the red button on the CareIP, in case of wall mounting the pull cord option can be used for activation. On activation, the carephone calls an alarm receiver to summon help.

There is a possibility to attach up to four different wired alarm triggers.

### **Different types of alarm transmitter**

Ten radio transmitters can be programmed with the CareIP, e.g. a radio trigger.

In the event of an alarm from a radio transmitter, the battery status of the radio transmitter is tested. The Elliot radio trigger has the option of automatic battery monitoring.

### **Alarm receivers**

Up to ten alarm receivers can be programmed. The alarm is identified using an alarm code which is programmed into the CareIP. An alarm placed to a private receiver is identified by a recorded voice ID. The person receiving the alarm can then speak directly with the person raising the alarm and take appropriate action.

### **Other functions**

Additional functions can be activated or changed by programming – refer to the service manual or contact your supplier for further information.

## Installation

### Installing the SIM card

- Unscrew the safety screw for the SIM card compartment cover.
- Remove the compartment cover.
- Install the SIM card according to picture below, note the markings on the card.
- Mount the cover and the safety screw for the SIM card compartment.

NOTE! 1234 or 0000 are not appropriate PIN codes for the SIM card, if any of these PIN codes are used, the code should be changed before installation. Change of PIN code is preferably done using a mobile phone.

### GSM signal strength check


Signal strength shall be checked with help of mobile phone with the same network operator as the SIM card installed in CareIP Mobile to find the best possible location of the unit.


## Connection

The CareIP Mobile is easy to install, if broadband/network connection shall be used, connection is made with the enclosed network cable.

Normally the CareIP is set to use DHCP, which means that the unit automatically will be assigned an IP address when connected to a network.

The CareIP should always be connected to the network as near incoming broadband as possible.

- Connect the power supply (1) to the wall socket and the power supply lead to the DC input socket .
- Connect the network cable (2) from the broadband/network to the Ethernet socket  on CareIP
- Turn the power switch  to the **I** position.
- Check that the “ON”  indicator lamp has a steady green light.
- The status light **+** shall go off within 30 sec. if GSM network is detected.
- Carry out a test alarm from the CareIP.

For further advice regarding installation, contact your supplier.

Installation must be performed only by authorized persons in accordance with the instructions in the service manual.

## Programming

Contacts, alarm codes and settings are programmed using a computer and the CareTech CS2 software or via an automatic download procedure from icare online.

In some cases, the CareIP can be delivered pre-programmed by the supplier.

For help with programming, please refer to the service manual or contact your supplier.

### Programming radio triggers

Up to ten radio triggers or other radio-alarm transmitters can be programmed to CareIP Mobile.

Programming radio transmitters:

- Hold down the green button for approximately four seconds until the green “On”  light starts to flash.
- Press the radio trigger. Once the CareIP receives the radio signal, you will hear a distinct tone signal.
- Confirm by pressing the red button on the carephone.

### Check radio coverage

This should always be done on installation.

- Hold down the green button for approximately four seconds until the green “On”  light starts to flash.
- Press the radio trigger. The built-in speaker will sound a tone signal on successful reception. Repeat this in all areas where the alarm needs to operate.

### Functional test




Always check the function after installation and programming. Test alarms shall be done regularly and in a way that secures the whole alarm chain from radio transmitter to alarm receiver.

## Radio trigger Elliot

Radio trigger Elliot uses two-way radio communication, this means that when CareIP Mobile receives an alarm from Elliot it will send an acknowledgement of the receipt. This is shown on Elliot by the indicator light that changes colour from red to green light. If Elliot doesn't get any acknowledgement the alarm will be repeated.

Elliot is equipped with two buttons that have programmable function

Accessories included: wristband, wristband attachment, neckband, pocket clip, "easy-press" accessory and cover. The accessories are attached to the alarm button with a simple click-in action.

## General information

- If Broadband/network connection is used the unit should be connected as close as possible to the incoming broadband.
- Verification of GSM coverage should always be done as a part of installation.
- In order for an alarm to be sent, at least one alarm receiver/response centre needs to be programmed.
- The carephone must not come into contact with water or any other liquid.
- Used batteries should be disposed of in an environmentally-responsible manner.

## Maintenance

### Cleaning

The carephone and radio trigger should be cleaned using a slightly damp cloth. Do not use solvents when cleaning.

### Replacing the CareIP battery

Battery should be replaced a.s.a.p after battery alarm or within three years.

### Replacing the radio trigger battery

To replace the radio trigger battery, please refer to the service manual or contact your supplier.

Battery replacement may be performed only by authorized persons in accordance with the instructions in service manual

## Technical data

### Carephone

Dimensions:  
Power supply:  
Backup battery:  
Backup time:  
Inputs:  
Output:  
Communication:  
Radio frequency:  
Radio receiver class:  
Number of radio transmitters:

### CareIP Mobile

200 x 175 x 35 mm (LxWxH)  
7,5 VDC power supply, (PoE optional)  
4,8V NiMH  
Up to two days  
Four inputs  
Three outputs, one Relay output, max 0,7A  
IP/(SIP), GSM  
869.2125 MHz  
1\*  
Up to 9

### GSM

Communication: GSM quad band 850/900/1800/1900 MHz  
GPRS class: 10  
SIM card interface: 1,8V and 3V

### Portable alarm button

Dimensions:  
Battery:  
Battery life:  
Water resistance:

### Elliot

45 x 30 x 13 mm (LxWxH)  
3 V lithium battery CR2430  
Up to five years  
Complies with IP67

### Equipment class:

Class 1 radio equipment

## Environment

This product is intended for indoor use in a normal residential environment.

Temperature: Operating temperature +5°C to +35°C  
Humidity: 0% to 75% relative humidity (non-condensed)  
Environmental class: 1 (EG-I)

\*Refers to radio module "C1", type 9340.

Radio receiver parameters: Max. usable sensitivity: -113dBm; Average usable sensitivity: ver. 28,7dBuV/m, hor.27,0dBuV/m; Adjacent channel selectivity: -26,2dBm/-28,2dBm; Intermodulation rejection: -36dBm/-37dBm; Spurious response rejection: -30dBm; Blocking or desensitisation: >-16dBm

CareIP can also be equipped with radio module type 9320, which allows maximum security and includes receipt of transmission, automatic transmission repeat, power control and frequency control. See order information for indicating which radio module the unit is equipped with.

# C€ 0402

Hereby, CareTech AB, Kalix Sweden declare that this equipment is in compliance with the essential requirements and other provisions of R&TTE Directive 1999/5/EC.

The comprehensive declaration of conformity is available at the address: CareTech AB, Box 10030, SE-952 27 Kalix, Sweden.

## Environmental information

This product complies with the requirements of the EU directives, 2002/96/EC (WEEE) and 2006/66/EG.

These directives regulate the product liability for electrical and battery recycling with the purpose of increasing recycling and minimising waste. The unit is marked with the "crossed out wheeled bin" logo, which indicates that it shall be handed in for recycling.

The product can be returned free of charge to a recycling station that is connected, directly or via a recycling system, to Care Tech or to your distributor. For detailed instructions, please check with your distributor or visit our website, [www.caretech.se](http://www.caretech.se)

Note! The WEEE information and recycling instructions applies to European Union member states only. For other countries please check local legislation or contact your distributor.

Manufactured in accordance with the EU directive, 2002/95/EC (RoHS)

The materials used in the radio trigger wristband and neckband meet the textile safety requirements of Oeko-Tex standard 100.

The information contained here is subject to change without notice. The only warranties for CareTech products and services are set forth in the express warranty statements accompanying such products and services. Nothing herein should be construed as constituting an additional warranty. CareTech shall not be liable for technical or editorial errors or omissions contained herein.

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